



No
disposable
electrodes
or cables
needed!

custo watch The new Holter Monitor



custo watch:

The result of consistent development work

The Holter is recorded with a lightweight chest belt, transmitted wirelessly to custo watch and is saved there. It comes without disposable electrodes or cables but with a high degree of patient comfort and an attractive design.

The new custo watch provides noticeable and measurable benefits for both the patient and the doctor.

Reliable

Powerful

Ergonomic

Top ECG quality

Modular

Precise



Chronometer

Date, Time



Status

Battery display, connection quality



Info

ECG view, heart rate



custo watch:

The advantages of a sophisticated system



custo belt 3

- Without adhesive electrodes
- Novel belt technology for 3-channel ECG recording
- Reusable
- Machine washable
- High wearing comfort
- With certificate of biocompatibility

If it is not possible to use the electrode belt, the ECG transmitter custo guard 3 can also be applied with an ECG cable and adhesive electrodes.



custo watch

- Signal testing of ECG channels
- Adjustable patient display
- Exchangeable bracelets
- Extremely lightweight
- Recording time: up to 72 h



custo guard 3

- 3-channel ECG transmitter
- Top ECG quality
- Extremely lightweight
- Energy-saving wireless transmission of the ECG signal to custo watch



custo docking station

- Charging of custo guard and custo watch
- Fast data transfer of Holter data to the PC
- Compatible with Holter evaluation software





custo holter ABDM:

Synchronous recording and view of Holter and ABPM


ABPM recorder custo screen 400

Beside the functions of custo screen 300, it additionally offers:

- Continuous saving of 3-channel ECG data from custo guard 3
- SD memory card for fast data download
- ABPM recording also possible without Holter

Software Holter and ABPM

- Combined view




Holter transmitter custo guard 3

- 3-channel ECG transmitter
- Top ECG quality
- Extremely lightweight
- Energy-saving wireless transmission of ECG signal to the ABPM recorder custo screen 400

3-channel electrode belt custo belt 3

- Holter without adhesive electrodes
- Novel belt technology for 3-channel ECG recording
- Reusable, machine washable
- High wearing comfort
- With certificate of biocompatibility



custo med – your local dealer

Contact:
 custo med GmbH
 Maria-Merian-Str. 6
 85521 Ottobrunn
 Germany

Phone: +49 89 71098 00
 Fax: +49 89 71098 10
 internet: www.customed.de
 e-mail: info@customed.de